

Index 2 Pulse Oximeter Simulator

Biomedical

Technical Data



The Index 2 is the most versatile optical simulator for oximeters on the market today. This lightweight, portable tool includes preloaded manufacturers' R-curves and the ability to define other "makes" for most pulse oximeters.

Motion presets, player mode, transmission level control (TLC), and computer commands boost testing ability. The Index 2 can also be configured to include an optional electrical simulation feature with probe test. Optical and electrical simulations allow technicians to isolate problems quickly. The probe test identifies defective probes with quantitative test results..

Key Features

- Portable
- 10 preloaded manufacturers' R-curves
- User-definable "makes" for most other manufacturers
- New R-curves for Masimo, Nonin and Philips Medical Systems (formerly Agilent / H-P) oximeters
- Six downloadable R-curve spaces available
- Simultaneous simulation of motion and arterial-oxygen levels
- Arterial wave-amplitude scale, calibrated in units of perfusion
- Tap/shiver motion simulations to explore the impact of motion
- RS232 port for computer control
- Physiological finger for complete SpO₂ tests
- Electrical simulations with probe testing (optional)

Technical Specifications

Range: 35 % to 100 % Resolution: 1 %

Accuracy: 100 % to 75 %: \pm 1 % \pm accuracy of the pulse oximeter under test; 74 % to 50 %: \pm 3 % \pm the accuracy of the pulse oximeter

under test; < 50 % unspecified Repeatability: ± 1 standard deviation

Rate

Range: 30 BPM to 250 BPM Resolution: 1 BPM Accuracy: 1 % ± 1 BPM

Pulse Amplitude

Range: 0 % to 100 % of nominal pleth amplitude

Resolution: 1 %

Pulse amplitude is 20 % of maximum

pass-through brightness

Probe Test

Continuity/Resistance

Test Matrix: Measures all combinations of possible interconnections in an XX point matrix

Range: 250 Ω to 150 k Ω Accuracy: \pm 5 % of reading

LED/Detector Voltage Test

Test Format: Measures the voltage drop across Red LED, infrared LED, and the photo detector when the internally generated test signal is

Test Signal: Constant current source @ 1 mA

Open Circuit: 2.5 V max

Measurement/Display Range: 0 V to 4 V Accuracy: ± 5 % of reading, 0.4 V to 4 V

Dynamic Test

Test Format: Photodetector/diode response to both the red and infrared light generated by the probe when pulsed by an internal test signal Test Signal: Pulsed between the two LEDs;

constant current level @ 1 mA

Test Results: Nominal range of 0 to 2000

About Fluke Biomedical
Fluke Biomedical is the world's leading manufacturer
of quality biomedical test and simulation products. In
addition, Fluke Biomedical provides the latest medical
imaging and oncology quality-assurance solutions for
regulatory compliance.

regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment
As a medical device manufacturer, we recognize and
follow certain quality standards and certifications when
developing our products. We are ISO 9001 certified and
our products are:

• FDA Compliant
• CE Certified, where required
• NIST Traceable and Calibrated
• UL, CSA, ETL Certified, where required

Checksum

Sum of all locations in the program chip; for service use only

General Information

Display: 2-line x 24-character super twist LCD

Battery Life: At least 4 hours of

continuous use

Dimensions: 18 in L x 16 in W x 9 in H (45.7 cm L x 40.6 cm W x 22.9 cm H)

Weight: 10 lb (4.5 kg)



Biomedical

Ordering Information

Model

2250232: Index 2-USA 120V 2395290: Index 2-SHK250V

Standard Accessories

2226196: Operator's manual

2213080: Battery charger (USA 120 Vac to

12 Vdc)

2213102: Battery charter (Euro/Schuko 220

Vac to 12 Vdc)

Optional Accessories

2248568: Multipurpose hard-sided

watertight carrying case

2204282: Soft vinyl carrying case

2200102: Interface cable, medTester to Index 2 (RS232; Female DB25 to

female DB9)

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical

PO Box 9090, Everett, WA 98206-9090 USA

Fluke Biomedical Europe AS

Vegamot 8, N-7048 Trondheim, Norway

For more information, contact us:

In the U.S.A. (800) 648-7952 or Fax (425) 446-5629

In Europe/M-East/Africa +47 73954700 or Fax +47 73954701

From other countries +1 (425) 347-6100 or Fax +1 (425) 446-5629

Email: sales@flukebiomedical.com

Web access: http://www.flukebiomedical.com

©2006 Fluke Biomedical. All rights reserved. Specifications subject to change without notice. Printed in U.S.A. LITO50